




Disposable Circular Haemorrhoid Stapler



CE 0473



Visible indicator makes it easier for the operator to keep the staple height range between 0.8 – 1.5mm



Unique hardened cutting ring provides the precise cutting of the tissue



Safe and solid closure of anvil and cartridge

Surgical Accessories



P02 anorectal introducer



P03 anorectal dilator



P04 anorectal transfixion tool



suture hook

Reference Code	Number of Staples	Diameter (mm)	Height (mm)
WM-CSH-32	30	32	4.2
WM-CSH-34	32	34	4.2