

DISPOSABLE HEALTHCARE PRODUCTS



## Urology

WORKING IN PARTNERSHIP WITH  
HEALTHCARE PROFESSIONALS

**Pennine**<sup>®</sup>  
H E A L T H C A R E



<b>Contents</b>		<b>Page</b>
<b>Catheters</b>	Nelaton Male Catheter	3
	Female Catheter	3
	Tiemann Catheter	4
<b>Pre-lubricated Catheters</b>	PreLube Catheters	5
	PreLube Bag	5
<b>Self Catheterisation System</b>	Pennine Intro-Uri-Cath	6
<b>Urine Collection</b>	Mid-Stream Urine Specimen Collector	7



## Nelaton (male) Catheters

For intermittent bladder drainage

### Nelaton (Male) Catheters

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- 100 pouches per inner polythene bag
- 38cm length
- Soft rounded tip for atraumatic insertion
- Two side eyes for maximum drainage
- Supplied single wrapped (SW) as standard

Single Wrapped Product Code	Qty Per Box	Double Wrapped Product Code	Qty Per Box	CH/Diam
NC-1206/SW	500	NC-1206/DW	400	6ch Ø 2.00mm
NC-1208/SW	500	NC-1208/DW	400	8ch Ø 2.66mm
NC-1210/SW	500	NC-1210/DW	400	10ch Ø 3.33mm
NC-1212/SW	500	NC-1212/DW	400	12ch Ø 4.00mm
NC-1214/SW	500	NC-1214/DW	400	14ch Ø 4.66mm
NC-1216/SW	500	NC-1216/DW	400	16ch Ø 5.33mm
NC-1218/SW	500	NC-1218/DW	400	18ch Ø 6.00mm
NC-1220/SW	400	NC-1220/DW	300	20ch Ø 6.66mm
NC-1222/SW	400	NC-1222/DW	300	22ch Ø 7.33mm
NC-1224/SW	400	NC-1224/DW	300	24ch Ø 8.00mm

### NHS supply chain codes

- Double wrapped (DW) to provide aseptic non-touch technique
- Ward pack configuration

Product Code	NHS NPC Code	Qty Per Box	CH/Diam
NC-1210/DWTRE	FST627	25	10ch Ø 3.33mm
NC-1212/DWTRE	FST628	25	12ch Ø 4.00mm
NC-1214/DWTRE	FST629	25	14ch Ø 4.66mm

## Female Catheters

For intermittent bladder drainage

### Female Catheters

- 23cm length
- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- 100 pouches per inner polythene bag
- Soft rounded tip for atraumatic insertion
- Two side eyes for maximum drainage

Product Code	Qty Per Box	CH/Diam
FC-1406	1000	6ch Ø 2.00mm
FC-1408	1000	8ch Ø 2.66mm
FC-1410	1000	10ch Ø 3.33mm
FC-1412	1000	12ch Ø 4.00mm
FC-1414	1000	14ch Ø 4.66mm
FC-1416	1000	16ch Ø 5.33mm

### NHS supply chain codes

- Ward pack configuration

Product Code	NHS NPC Code	Qty Per Box	CH/Diam
FC-1406TRE	FST626	25	6ch Ø 2.00mm
FC-1408TRE	FST022	25	8ch Ø 2.66mm
FC-1410TRE	FST023	25	10ch Ø 3.33mm
FC-1412TRE	FST024	25	12ch Ø 4.00mm

Please note: All dimensions are nominal

Tel: 44 (0) 1332 794880 Email: [marketing@penninehealthcare.co.uk](mailto:marketing@penninehealthcare.co.uk) [www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)



## Tiemann Catheter

Traditional Nelaton with angled, pointed tip for directional stability

- Slightly hardened tip for ease of insertion through areas of obstruction
- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel packs, single wrapped
- Single use only
- 38cm length

Product Code	Qty Per Box	CH/Diam
NC-1708	100	8ch Ø 2.66mm
NC-1710	100	10ch Ø 3.33mm
NC-1712	100	12ch Ø 4.00mm
NC-1714	100	14ch Ø 4.66mm
NC-1716	100	16ch Ø 5.33mm
NC-1718	100	18ch Ø 6.00mm



## PreLube Catheters

The Pennine PreLube is a pre-lubricated catheter for Intermittent Urethral Self-Catheterisation. The PreLube is coated with Lubrajel® to reduce friction and minimise discomfort for the patient. The unique coating eradicates the need for water or any other gelling agent. The catheter is completely enclosed, making the PreLube very safe to use because there is no need to touch the catheter itself. The PreLube is ideal for teaching self-catheterisation.

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel packs, double wrapped
- Sealed pre-cut polythene wrapper
- Single use only
- Smooth closed end, two staggered eyes
- Hard PVC capsule containing Lubrajel® to offer total pre-lubrication of the Nelaton catheter and ensure quick and comfortable insertion
- Detachable packaging tip for aseptic introduction of the catheter into the urethra

Male PreLube (43cm)	Qty	Drug	CH/Diam
Product Code	Per Box	Tariff Code	
NC-1608	100	NC-1608/FP	8ch Ø 2.66mm
NC-1610	100	NC-1610/FP	10ch Ø 3.33mm
NC-1612	100	NC-1612/FP	12ch Ø 4.00mm
NC-1614	100	NC-1614/FP	14ch Ø 4.66mm
NC-1616	100	NC-1616/FP	16ch Ø 5.33mm

Female PreLube (23cm)	Qty	Drug	CH/Diam
Product Code	Per Box	Tariff Code	
FC-1608	100	FC-1608/FP	8ch Ø 2.66mm
FC-1610	100	FC-1610/FP	10ch Ø 3.33mm
FC-1612	100	FC-1612/FP	12ch Ø 4.00mm
FC-1614	100	FC-1614/FP	14ch Ø 4.66mm
FC-1616	100	FC-1616/FP	16ch Ø 5.33mm

- 1 litre collection bag
- 43cm Nelaton catheter
- Easily disposable

## PreLube Bag for PreLube Catheter

Product Code	Qty	Drug	CH/Diam
Code	Per Box	Tariff Code	
NC-1608/BG	100	NC-1608/BG/FP	8ch Ø 2.66mm
NC-1610/BG	100	NC-1610/BG/FP	10ch Ø 3.33mm
NC-1612/BG	100	NC-1612/BG/FP	12ch Ø 4.00mm
NC-1614/BG	100	NC-1614/BG/FP	14ch Ø 4.66mm
NC-1616/BG	100	NC-1616/BG/FP	16ch Ø 5.33mm

PreLube and PreLube Bags are available in the U.K. on Drug Tariff



## Pennine Intro-Uri-Cath

Sterile Field Urinary Catheter Introducer. A complete intermittent urethral catheterisation procedure kit. Extends the sterile field into the patient.

- Easy, quick, time saving, economical, safe evidence of reduced bacteriuria
- Sterile, individual peel pack with clear instructions for use
- Standard Nelaton catheter
- Silicone introducer with lubricant
- Graduated bag for dual purpose no-touch technique and collection of urine

Product Code	Qty Per Box	Drug Tariff Code	CH/Diam
NC-1508	20	NC-1508/FP	8ch Ø 2.66mm
NC-1510	20	NC-1510/FP	10ch Ø 3.33m
NC-1512	20	NC-1512/FP	12ch Ø 4.00mm
NC-1514	20	NC-1514/FP	14ch Ø 4.66mm

### Indications

#### Hospitals

- Where intermittent catheterisation is preferable to use of indwelling catheters
- Post trauma and post surgical bladder relief
- Prior to lower abdominal surgery
- Acute retention of urine and all situations where there is inadequate bladder drainage
- Spinal injuries, Spina Bifida where there is risk of U.T.I.
- Teaching self catheterisation
- Elderly and other patient groups where there is a susceptibility to nosocomial U.T.I.

#### Maternity

- Pre and post delivery bladder drainage to reduce chances of bacteriuria
- Measurement of residual urine

#### Rehabilitation/Community

- For teaching self catheterisation
- For many patient groups where long term intermittent catheterisation is associated with a need for a complete, self contained system and which reduces the risk of induced U.T.I
- For general practitioner and other community medical personnel for emergency use, e.g. urinary retention





## Mid-Stream Urine Specimen Collector

For aseptic collection of midstream urine

- Easy to use, 30ml container
- Patient details label
- Patient instructions

Product Code	Type	Trap Size	Qty Per Box
MSU-1000	Without cleansing wipe	30ml	50
MSU-1001	With cleansing wipes	30ml	50

Pennine Healthcare is one of the UK's largest privately owned manufacturers of single-use sterile medical devices, with a customer base, not only in the UK but in over 50 countries worldwide.

Pennine has an extensive range of disposables for use in urology, anaesthetics, respiratory, gastroenterology and surgical suction, as well as being a leading supplier of Custom Procedure packs.

Pennine's vision is to be a "world class, innovative manufacturer, supplying the global healthcare market with outstanding commitment to its people, customers and suppliers".

Pennine is enormously proud of its achievements but knows that it is only successful because of its people and real values, which are at the heart of everything they do

To be progressive  
To be responsive  
To be passionate  
To have integrity

With a culture and history of innovation and engineering expertise, Pennine is still a competitive supplier of disposables on a worldwide basis and today embraces the skills and techniques of 'Lean Sigma' manufacture, ensuring that efficiency and waste reduction are targeted and that Pennine staff are highly trained, motivated and entirely focussed on the customer.

Pennine is consistently pursuing new opportunities and embracing innovative solutions, working with clinicians and healthcare professionals to develop new products for advanced patient care.

For further information, please call our sales team on +44 (0) 1332 794880 or email us at: [marketing@penninehealthcare.co.uk](mailto:marketing@penninehealthcare.co.uk). Alternatively, you can view all our product ranges at our website [www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)



Gastroenterology



Microsuction Cannulae



Custom Procedure Packs



Respiratory Products



Surgical Suction



Airways Management

**Pennine**<sup>®</sup>  
H E A L T H C A R E

**Pennine Healthcare**

City Gate  
London Road  
Derby DE24 8WY  
United Kingdom  
Telephone: +44 (0) 1332 794880  
Fax: +44 (0) 1332 794890

Email: [info@penninehealthcare.co.uk](mailto:info@penninehealthcare.co.uk)  
Web: [www.penninehealthcare.co.uk](http://www.penninehealthcare.co.uk)

Some of the products in this catalogue are subject to a minimum order quantity.  
Please ask our sales department for further details.

Ivor Shaw Ltd trading as Pennine Healthcare