



DISPOSABLE HEALTHCARE PRODUCTS



Gastroenterology

WORKING IN PARTNERSHIP WITH
HEALTHCARE PROFESSIONALS

Pennine[®]
HEALTHCARE



Contents		Page
Catheters	Rectal Catheters	3
	Rectal Tubes	3
	Stomach Tubes	3
	Oesophageal Tubes	3
Ryles and Levins Tubes	Ryles Tubes	4
	Ryles Tubes (Metal Tipped)	4
	Levins/ Duodenal Tube	4
Bile Bags Catheter Spigots	Bile Bags	5
	Catheter Spigots	5
Infant Feeding Tubes	Standard	6
	Graduated with X-Ray line	6



Rectal Catheters 40cm

For atraumatic drainage and irrigation

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- Large side eyes and closed tip

Rectal Catheters



Large side eyes/closed tip

Code No.	Qty Per Box	CH/Diam
RC-1818	500	18ch Ø 6.00mm
RC-1820	400	20ch Ø 6.66mm
RC-1822	400	22ch Ø 7.33mm
RC-1824	400	24ch Ø 8.00mm
RC-1826	300	26ch Ø 8.66mm
RC-1828	300	28ch Ø 9.33mm
RC-1830	300	30ch Ø 10.00mm

Rectal Tubes 78cm

As above but with a wide open end for administration of enema solutions, drainage and irrigation

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches

Rectal Tubes



Wide open end

Code No.	Qty Per Box	CH/Diam
RT-2618	100	18ch Ø 6.00mm
RT-2620	50	20ch Ø 6.66mm
RT-2622	50	22ch Ø 7.33mm
RT-2624	50	24ch Ø 8.00mm
RT-2626	50	26ch Ø 8.66mm
RT-2628	50	28ch Ø 9.33mm
RT-2630	50	30ch Ø 10.00mm

Stomach/Oesophageal Tubes 78cm

Mainly for gastric lavage-washout-or other gastric surgical procedures. With or without X-ray detectable line.

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- 2 side eyes, closed end
- Optional with X-ray opaque line, XST type

Stomach Tubes



Two side eyes/closed tip

Stomach Tubes Code No.	Stomach Tubes x-ray line Code No.	Qty Per Box	CH/Diam
ST-2408	XST-2408	100	8ch Ø 2.66mm
ST-2410	XST-2410	100	10ch Ø 3.33mm
ST-2412	XST-2412	100	12ch Ø 4.00mm
ST-2414	XST-2414	100	14ch Ø 4.66mm
ST-2416	XST-2416	100	16ch Ø 5.33mm
ST-2418	XST-2418	100	18ch Ø 6.00mm
ST-2420	XST-2420	50	20ch Ø 6.66mm
ST-2422	XST-2422	50	22ch Ø 7.33mm
ST-2424	XST-2424	50	24ch Ø 8.00mm
ST-2426	XST-2426	50	26ch Ø 8.66mm
ST-2428	XST-2428	50	28ch Ø 9.33mm
ST-2430	XST-2430	50	30ch Ø 10.00mm



Ryles Tubes 122cm

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- Printed indications at 40-50-60 and 70cm from tip
- With X-Ray opaque tip, 4 side eyes, closed end
- With or without removable adaptor and integral cap (Code number /L)
- Optional of X-Ray opaque tip, X-Ray opaque tip + Luer available

Ryles Tubes



With x-ray tip Code No.	With Luer & x-ray tip Code No.	With x-ray line/ Tip & Luer Code No.	Qty Per Box	CH/Diam
RT-2008	RT-2008/L	RT-2308	100	8ch Ø 2.66mm
RT-2010	RT-2010/L	RT-2310	100	10ch Ø 3.33mm
RT-2012	RT-2012/L	RT-2312	100	12ch Ø 4.00mm
RT-2014	RT-2014/L	RT-2314	100	14ch Ø 4.66mm
RT-2016	RT-2016/L	RT-2316	100	16ch Ø 5.33mm
RT-2018	RT-2018/L	RT-2318	100	18ch Ø 6.00mm
RT-2020		RT-2320	50	20ch Ø 6.66mm
RT-2022		RT-2322	50	22ch Ø 7.33mm

Metal Tipped Ryles Tubes 122cm

As above but with stainless Steel weighted tip to aid intubation and positioning in patient

- X-ray opaque line
- Printed indications at 40-50-60 and 70cm from tip
- Loose Luer adaptor included in the pack

Metal Tip Ryles Tubes



Code No.	Qty Per Box	CH/Diam
RT-2108	100	8ch Ø 2.66mm
RT-2110	100	10ch Ø 3.33mm
RT-2112	100	12ch Ø 4.00mm
RT-2114	100	14ch Ø 4.66mm
RT-2116	100	16ch Ø 5.33mm
RT-2118	100	18ch Ø 6.00mm

Levins/Duodenal Tube 122cm

Duodenal tube for active aspiration, feeding, injection, lavage or other gastric surgical procedures. With or without X-ray detectable line

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- Printed indications at 40-50-60 and 70cm from tip
- 4 side eyes to reduce risk of blockage
- Terminal eye to ensure maximum dispersion
- Optional closed end
- Optional with X-ray opaque line, XLT type

Levins Tubes, no X-Ray line



Levins Tubes, X-Ray line



Open end Code No.	Open end & x-ray line Code No.	Closed end Code No.	Closed end & x-ray line Code No.	Qty Per Box	CH/Diam
LT-2206	XLT-2206	LT-2406	XLT-2406	100	6ch Ø 2.00mm
LT-2208	XLT-2208	LT-2408	XLT-2408	100	8ch Ø 2.66mm
LT-2210	XLT-2210	LT-2410	XLT-2410	100	10ch Ø 3.33mm
LT-2212	XLT-2212	LT-2412	XLT-2412	100	12ch Ø 4.00mm
LT-2214	XLT-2214	LT-2414	XLT-2414	100	14ch Ø 4.66mm
LT-2216	XLT-2216	LT-2416	XLT-2416	100	16ch Ø 5.33mm
LT-2218	XLT-2218	LT-2418	XLT-2418	100	18ch Ø 6.00mm
LT-2220	XLT-2220	LT-2420	XLT-2420	50	20ch Ø 6.66mm

Please note: All angles and lengths are nominal

Tel: 44 (0) 1332 794880 Email: marketing@penninehealthcare.co.uk www.penninehealthcare.co.uk

LATEX
FREE



Bile Bags

For collection of Bile and other gastric secretions. Easy connection to a Pennine nasogastric tube (Ryles tube or Levins tube) or with Kehr's tube after Bile duct surgery

- Choice of two bottom outlets; capped or stop cock for cleaner emptying procedure
- Available with latex straps for bed or patient attachment

Code No.	Top Connector	Bottom Outlet	Qty
DB3B with straps	Kehr's tube or female luer connection	Capped Outlet	100
DB3B/R with straps	Ryles/Levins tube connection	Capped Outlet	100
BB350T/Y	Kehr's tube or female luer connection	Stop cock	100
BB350R/Y	Ryles/Levins tube connection	Stop cock	100


Catheter Spigots

A range of sterile tubing closures for catheters and tubes

- Sterile
- Individual peel pouches

Code No.	Size	Length	Max. diam.	Qty
CS-9001	Small	56mm	6mm	500
CS-9002	Medium	68mm	9mm	500
CS-9003	Large	78mm	11mm	500





Infant Feeding Tubes

40cm, 48cm, 105cm - For enteral feeding of small babies and infants

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel pouches
- 100 tubes per polythene bag
- Closed end, two side eyes
- Printed graduations at 10 and 20cm from tip
- Standard, plain PVC

40cm Code No.	48cm Code No.	Qty Per Box	105cm Code No.	Qty Per Box	CH/Diam
FT-1605/40	FT-1605/48	500	FT-1605/105	100	5ch Ø 1.66mm
FT-1606/40	FT-1606/48	500	FT-1606/105	100	6ch Ø 2.00mm
FT-1608/40	FT-1608/48	500	FT-1608/105	100	8ch Ø 2.66mm
FT-1610/40	FT-1610/48	500	FT-1610/105	100	10ch Ø 3.33mm

Graduated Infant Feeding Tubes with X-ray lines

40cm, 48cm, 105cm

- Medical grade PVC (Shore Hardness A80)
- Sterile, individual peel packs
- 40cm, 48cm or 105cm long, closed end - 2 side eyes
- 40cm, 48cm length printed indications at every cm between 5cm and 25cm from tip
- 105cm length printed indications at every cm between 5cm and 70cm from tip
- X-ray opaque line

40cm Code No.	48cm Code No.	Qty Per Box	105cm Code No.	Qty Per Box	CH/Diam
FT-1905/40	FT-1905/48	500	FT-1905/105	100	5ch Ø 1.66mm
FT-1906/40	FT-1906/48	500	FT-1906/105	100	6ch Ø 2.00mm
FT-1908/40	FT-1908/48	500	FT-1908/105	100	8ch Ø 2.66mm
FT-1910/40	FT-1910/48	500	FT-1910/105	100	10ch Ø 3.33mm



Standard, plain PVC



Graduated X-ray tube

Please note: All angles and lengths are nominal

Tel: 44 (0) 1332 794880 Email: marketing@penninehealthcare.co.uk www.penninehealthcare.co.uk

LATEX
FREE 7

Pennine Healthcare is one of the UK's largest privately owned manufacturers of single-use sterile medical devices, with a customer base, not only in the UK but in over 50 countries worldwide.

Pennine has an extensive range of disposables for use in urology, anaesthetics, respiratory, gastroenterology and surgical suction, as well as being a leading supplier of Custom Procedure packs.

Pennine's vision is to be a "world class, innovative manufacturer, supplying the global healthcare market with outstanding commitment to its people, customers and suppliers".

Pennine is enormously proud of its achievements but knows that it is only successful because of its people and real values, which are at the heart of everything they do

To be progressive
To be responsive
To be passionate
To have integrity

With a culture and history of innovation and engineering expertise, Pennine is still a competitive supplier of disposables on a worldwide basis and today embraces the skills and techniques of 'Lean Sigma' manufacture, ensuring that efficiency and waste reduction are targeted and that Pennine staff are highly trained, motivated and entirely focussed on the customer.

Pennine is consistently pursuing new opportunities and embracing innovative solutions, working with clinicians and healthcare professionals to develop new products for advanced patient care.

For further information, please call our sales team on +44 (0) 1332 794880 or email us at: marketing@penninehealthcare.co.uk. Alternatively, you can view all our product ranges at our website www.penninehealthcare.co.uk



Microsuction Cannulae



Urology



Custom Procedure Packs



Respiratory Products



Surgical Suction



Airways Management

Pennine[®]
H E A L T H C A R E

Pennine Healthcare

City Gate
London Road
Derby DE24 8WY
United Kingdom
Telephone: +44 (0) 1332 794880
Fax: +44 (0) 1332 794890

Email: info@penninehealthcare.co.uk
Web: www.penninehealthcare.co.uk

Some of the products in this catalogue are subject to a minimum order quantity.
Please ask our sales department for further details.

Ivor Shaw Ltd trading as Pennine Healthcare