

DISPOSABLE HEALTHCARE PRODUCTS



Airway Management

WORKING IN PARTNERSHIP WITH
HEALTHCARE PROFESSIONALS

Pennine[®]
HEALTHCARE



Contents

Page

Catheters

Suction Catheters	3
Prestrol Suction Catheters	4
Finger Tip Vacuum Control Catheters	4
Oppo-Cath Suction Catheters	5
Prestrol Oppo-Cath Suction Catheters	5
Coude Tip Oppo-Cath Suction Catheters	6
Paediatric Suction Catheters	6
Clear Graduated Oppo-Cath Suction Catheters	7

Mucus Extractors

Mucus Specimen Traps	8
Mucus Extractors	8
Economy Mucus Extractor	9
Safe-Ex Mucus Extractor	9

Oxygen Therapy

Guedel Airways	10
Oxygen Connecting Tubing	10
Link Oxygen Connecting Tubing	10
Oxygen Catheters	11
Twin Nasal Oxygen Cannulae	11



Suction Catheters

48cm/53cm/60cm

For oral, tracheal and endobronchial suction procedures. To be used by connection to the suction source with Pennine suction connecting tubing, MF, FFM, or FFM/VC type.

- PVC Medical grade
- Sterile individual peel pouches
- Colour coded funnel end
- Available in softer grade PVC=/S
- Standard configuration=
terminal eye, 1 side eye
- Double offset eye option=/DE
- Terminal eye only option=/TO



48cm	48cm Soft	48cm Double Eye	48cm Terminal Eye Only	Qty. Per Box	CH/Dia.
SC-1005	-	-	-	500	5ch Ø 1.66mm
SC-1006	-	SC-1006/DE	SC-1006/TO	500	6ch Ø 2.00mm
SC-1007	-	-	-	500	7ch Ø 2.33mm
SC-1008	-	SC-1008/DE	SC-1008/TO	500	8ch Ø 2.66mm
SC-1010	SC-1010/S	SC-1010/DE	SC-1010/TO	500	10ch Ø 3.33mm
SC-1012	SC-1012/S	SC-1012/DE	SC-1012/TO	500	12ch Ø 4.00mm
SC-1014	SC-1014/S	SC-1014/DE	SC-1014/TO	500	14ch Ø 4.66mm
SC-1016	SC-1016/S	SC-1016/DE	SC-1016/TO	500	16ch Ø 5.33mm
SC-1018	SC-1018/S	SC-1018/DE	SC-1018/TO	500	18ch Ø 6.00mm
SC-1020	-	-	-	400	20ch Ø 6.66mm

53cm	53cm Double Eye	53cm Terminal Eye Only	Qty. Per Box	CH/Dia.
SC-1006/53	-	-	-	6ch Ø 2.00mm
SC-1008/53	-	SC-1008/53/TO	500	8ch Ø 2.66mm
SC-1010/53	SC-1010/53/DE	SC-1010/53/TO	500	10ch Ø 3.33mm
SC-1012/53	SC-1012/53/DE	SC-1012/53/TO	500	12ch Ø 4.00mm
SC-1014/53	SC-1014/53/DE	SC-1014/53/TO	500	14ch Ø 4.66mm
SC-1016/53	SC-1016/53/DE	SC-1016/53/TO	500	16ch Ø 5.33mm
SC-1018/53	SC-1018/53/DE	-	500	18ch Ø 6.00mm

60cm	60cm Double Eye	Qty. Per Box	CH/Dia.
SC-1010/60	SC-1010/DE/60	100	10ch Ø 3.33mm
SC-1012/60	SC-1012/DE/60	100	12ch Ø 4.00mm
SC-1014/60	SC-1014/DE/60	100	14ch Ø 4.66mm
SC-1016/60	SC-1016/DE/60	100	16ch Ø 5.33mm
SC-1018/60	-	100	18ch Ø 6.00mm



Prestrol Suction Catheters

48cm/53cm/60cm

Vacuum controlled Suction Catheter for oral, tracheal and endobronchial suction procedures. To be used by connection to the suction source with Pennine suction connecting tubing, FF type.

- PVC Medical Grade
- Sterile individual peel pouches
- Colour coded Prestrol end enabling vacuum control
- Standard configuration = terminal eye/1 side eye
- 100 per polythene bag



48cm	Qty. Per Box	60cm	Qty. Per Box	CH/Dia.
SC-1105	400	-	-	5ch Ø 1.66mm
SC-1106	400	SC-1106/60	100	6ch Ø 2.00mm
SC-1108	400	SC-1108/60	100	8ch Ø 2.66mm
SC-1110	400	SC-1110/60	100	10ch Ø 3.33mm
SC-1112	400	SC-1112/60	100	12ch Ø 4.00mm
SC-1114	400	SC-1114/60	100	14ch Ø 4.66mm
SC-1116	400	SC-1116/60	100	16ch Ø 5.33mm
SC-1118	400	-	-	18ch Ø 6.00mm

Fingertip Vacuum Control Catheters

- Double side eyes=DE
- With capped vacuum control connector



VC2 Control Length 53cm	Qty. Per Box	CH/Dia.
SC-2006/53/DE	100	6ch Ø 2.00mm
SC-2008/53/DE	100	8ch Ø 2.66mm
SC-2010/53/DE	100	10ch Ø 3.33mm
SC-2012/53/DE	100	12ch Ø 4.00mm
SC-2014/53/DE	100	14ch Ø 4.66mm
SC-2016/53/DE	100	16ch Ø 5.33mm
SC-2018/53/DE	100	18ch Ø 6.00mm



Oppo-Cath Suction Catheters

40cm/53cm

Suction catheters for oral, tracheal and endobronchial suction procedures. Two side eyes for increased patient safety and reduced tissue damage. To be used by connection to the suction source with Pennine suction connecting tubing, MF type.

- Terminal eye, two small opposed side eyes
- Funnel end
- Sterile, individual peel pouch
- 100 units per polythene bag



40cm	53cm	Qty. Per Box	53cm	Qty. Per Box	CH/Dia.
OSC-1306/40	OSC-1306/53	500	OSC-1306/53/1	100	6ch Ø 2.00mm
OSC-1308/40	OSC-1308/53	500	OSC-1308/53/1	100	8ch Ø 2.66mm
OSC-1310/40/H	OSC-1310/53/H	500	OSC-1310/53/H/1	100	10ch Ø 3.33mm
OSC-1312/40/H	OSC-1312/53/H	500	OSC-1312/53/H/1	100	12ch Ø 4.00mm
OSC-1314/40/H	OSC-1314/53/H	500	OSC-1314/53/H/1	100	14ch Ø 4.66mm
OSC-1316/40/H	OSC-1316/53/H	500	OSC-1316/53/H/1	100	16ch Ø 5.33mm
-	OSC-1318/53/H	500	OSC-1318/53/H/1	100	18ch Ø 6.00mm



Prestrol Oppo-Cath Suction Catheters

48cm/60cm

Vacuum control Suction Catheter for oral, tracheal and endobronchial suction procedures. Two side eyes for increased patient safety and reduced tissue damage. To be used by connection to the suction source with Pennine suction connecting tubing, FF type.

- PVC Medical grade
- Terminal eye, two small opposed side eyes
- Sterile, individual peel pouches
- 100 units per polythene bag
- With colour coded Prestrol Vacuum Control device



48cm	Qty. Per Box	60cm	Qty. Per Box	CH/Dia.
OSC-1505/48	400	-	100	5ch Ø 1.66mm
OSC-1506/48	400	-	100	6ch Ø 2.00mm
OSC-1508/48	400	OSC-1508/60	100	8ch Ø 2.66mm
OSC-1510/48	400	OSC-1510/60	100	10ch Ø 3.33mm
OSC-1512/48	400	OSC-1512/60	100	12ch Ø 4.00mm
OSC-1514/48	400	OSC-1514/60	100	14ch Ø 4.66mm
OSC-1516/48	400	OSC-1516/60	100	16ch Ø 5.33mm
OSC-1518/48	400	OSC-1518/60	100	18ch Ø 6.00mm





Coude Tip Oppo-Cath Suction Catheters

53cm/60cm

Identical to our Oppo-Cath Suction Catheters, but with an angled tip. For endobronchial suction procedures. Gives easier access to the right or left bronchus. To be used by connection to the suction source using Pennine suction connecting tubing, MF, FFM or FFM/VC type.

- Funnel end
- Sterile, individual peel pouch
- 100 units per polythene bag
- Terminal eye, two small opposed side eyes

53cm	Qty. Per Box	60cm	Qty. Per Box	CH/Dia.
OSC-1710/53/H	400	-	100	10ch Ø 3.33mm
OSC-1712/53/H	400	OSC-1712/60/H	100	12ch Ø 4.00mm
OSC-1714/53/H	400	OSC-1714/60/H	100	14ch Ø 4.66mm
OSC-1716/53/H	400	OSC-1716/60/H	100	16ch Ø 5.33mm



Paediatric Suction Catheters

25cm/30cm

For oral, tracheal and bronchial suctioning on babies and children.

- Medical grade PVC
- Sterile, individual peel pouches

Oppo-Cath Suction Catheters (OSC-13 range)

- 25cm
- Standard colour coded funnel ends
- 2 small opposed side eyes

Xray Graduated Suction Catheters (OSC-14 range)

- 30cm
- Xray line
- Graduation every 1cm
- Colour coded luer connection
- 2 small opposed side eyes

Prestrol Oppo-Cath (OSC-16 range)

- 30cm
- Colour coded vacuum control connector
- 2 small opposed side eyes

Oppo Cath 25cm	Qty. Per Box	Xray Graduated 30cm	Prestrol Oppo Cath 30cm	Qty. Per Box	CH/Dia.
OSC-1306/25	500	OSC-1406/30/G	OSC-1606/30	100	6ch Ø 2.00mm
OSC-1308/25	500	OSC-1408/30/G	OSC-1608/30	100	8ch Ø 2.66mm
OSC-1310/25/H	500	OSC-1410/30/G	OSC-1610/30	100	10ch Ø 3.33mm
OSC-1312/25/H	500	-	-	-	12ch Ø 4.00mm
OSC-1314/25/H	500	-	-	-	14ch Ø 4.66mm



Clear Graduated Oppo-Cath Suction Catheters

48cm/60cm

Suction catheters for oral, tracheal and endobronchial suction procedures.

For increased patient safety and reduced tissue damage. To be used by connecting to the suction source with Pennine suction connecting tubing, FF type.

- PVC Medical grade
- Clear Prestrol Vacuum control
- Graduations every 1cm
- Sterile, individual peel pouch
- 100 units per polythene bag
- Terminal eye, two small opposed side eyes

Clear Prestrol 48cm	Qty. Per Box	Clear Prestrol 60cm	Qty. Per Box	CH/Dia.
OSC-1205/48	400	-	-	5ch Ø 1.66mm
OSC-1206/48	400	-	-	6ch Ø 2.00mm
OSC-1207/48	400	-	-	7ch Ø 2.33mm
OSC-1208/48	400	OSC-1208/60	100	8ch Ø 2.66mm
OSC-1210/48	400	OSC-1210/60	100	10ch Ø 3.33mm
OSC-1212/48	400	OSC-1212/60	100	12ch Ø 4.00mm
OSC-1214/48	400	OSC-1214/60	100	14ch Ø 4.66mm
OSC-1216/48	400	OSC-1216/60	100	16ch Ø 5.33mm





Mucus Specimen Traps

To collect mucus specimen for micro-biological examination.

- Sterile, individual pack
- 25ml or 70ml trap
- Attachment for suction catheter and suction connecting tube
- Patient details label
- Spare cap, screw fit

Product Code	Capacity	Type	Qty. Per Box
MST-3000	25ml	With spare screw cap	50
MST-3070	25ml	With reduced length connection to suction catheter to ease movement of secretions to the trap. With spare screw cap & patient label	50
MST-4000	70ml	For use in bronchoscopy. With spare screw cap	50



MST-3000



MST-3070



MST-4000

Mucus Extractors

Mucus extractor designed to obtain mucus specimen for micro-biological examination with 33cm Catheter.

- 20ml Trap
- With colour coded cap for easy identification of catheter size
- With standard funnel or vacuum control connector for attachment to suction source

Without Vacuum Control Connector

Product Code	Qty. Per Box	Catheter Size
ME-8108	75	8ch Ø 2.66mm
ME-8110	75	10ch Ø 3.33mm
ME-8112	75	12ch Ø 4.00mm
ME-8114	75	14ch Ø 4.66mm



ME-8108

With Vacuum Control Connector

Product Code	Qty. Per Box	Catheter Size
ME-9108	75	8ch Ø 2.66mm
ME-9110	75	10ch Ø 3.33mm
ME-9112	75	12ch Ø 4.00mm
ME-9114	75	14ch Ø 4.66mm



ME-9108



Economy Mucus Extractor

Easy to use, familiar standard system.

Available with female fitting for mechanical suctioning.

- 25ml Trap
- 14ch, 33cm suction catheter with terminal eye only
- Soft PVC, sterile, straight peel pouch

Product Code	Trap Capacity	Qty. Per Box	Catheter Size
ME-6004	25ml	100	14ch Ø 4.66mm



Safe-Ex Mucus Extractors

Unique, patented mucus extraction using the 'Safe-Ex' block filter device. Provides maximum protection for user by irreversible blocking on contact with any aqueous fluids such as mucus. Microbiologically & clinically tested. Easy to use, familiar standard system using mouth suctioning procedures.

With sample collection 25ml trap with extra screw cap for transport to the laboratory.

- 25ml Trap
- 10ch or 14ch suction catheter
- Spare cap, screw fit

Product Code	Trap Capacity	Qty. Per Box	Catheter Size
ME-7100	25ml	50	10ch Ø 3.33mm
ME-7140	25ml	50	14ch Ø 4.66mm





Guedel Airways

Oropharyngeal airways for use during anaesthesia or emergency.


- Sterile, individual peel pouches or bulk, non-sterile
- Colour-coded bite-blocks for easy size identification
- Manufactured in linear low-density polyethylene

Product Code (Sterile)	Qty. Per Box	Product Code (Bulk)	Qty. Per Box	Length	Size
GA-92000	100	GA-92000/B	500	34mm	000
GA-9200	100	GA-9200/B	500	45mm	00
GA-9210	100	GA-9210/B	500	55mm	0
GA-9211	100	GA-9211/B	500	71mm	1
GA-9212	100	GA-9212/B	500	90mm	2
GA-9213	100	GA-9213/B	400	100mm	3
GA-9214	100	GA-9214/B	300	110mm	4
GA-9215	50	GA-9215/B	250	115mm	5


Oxygen Connecting Tubing

For connecting oxygen source to oxygen administration device.

- Clean, non-sterile
- Individual polythene packs
- Female funnel connectors to fit most oxygen therapy devices
- Also available as kink resistant Star Lumen tubing



Standard Tubing



Star Lumen Tubing

- Clean, non-sterile
- Individual polythene packs
- Link female, funnel ends can be easily cut to provide sure fit to any oxygen therapy device
- Also available in 30m rolls



Standard Tubing

Product Code	Qty. Per Box	Length
OC-7000	100	1.8m
OC-7009	100	2.7m

Star Lumen Tubing


Product Code	Qty. Per Box	Length
OC-9000	100	1.8m
OC-9009	100	2.4m

Link Oxygen Connecting Tubing (Bubble Tubing)

Product Code	Qty. Per Box	Length
OC-7311	100	1.8m
OC-7341	100	2.7m

Rolls in Dispenser Boxes

Product Code	Qty. Per Box	Length
OL-8400	1 Roll	30m



Oxygen Catheters

Catheters for direct nasal administration of oxygen.
Male fitting connects to oxygen source.

- Sterile, individual peel pouch
- 100 catheters per polythene bag
- Green PVC
- Male connectors
- Available with closed end, multiple side holes or with a sponge compress tip

With holes		
Product Code	Qty. Per Box	CH/Dia
OC-7010	500	10ch Ø 3.33mm
OC-7012	500	12ch Ø 4.00mm
OC-7014	500	14ch Ø 4.66mm

With Sponge Compress		
Product Code	Qty. Per Box	CH/Dia
OC-7110	500	10ch Ø 3.33mm
OC-7112	500	12ch Ø 4.00mm
OC-7114	500	14ch Ø 4.66mm

- Link female funnel connection
- 220cm length
- Multiple side eyes

Link Oxygen Catheter with Holes - 220cm		
Product Code	Qty. Per Box	CH/Dia
OC-7210	100	10ch Ø 3.33mm



Twin Nasal Oxygen Cannulae

- Clean, non-sterile
- Individual polythene bag
- Green PVC
- Twin nasal cannulae with nasal prongs and adjustable head piece
- 2m of connecting tubing with female funnel connections

Product Code	Qty. Per Box	Length
ON-7200	100	2.0m



Pennine Healthcare is one of the UK's largest privately owned manufacturers of single-use sterile medical devices, with a customer base, not only in the UK but in over 50 countries worldwide.

Pennine has an extensive range of disposables for use in urology, anaesthetics, respiratory, gastroenterology and surgical suction, as well as being a leading supplier of Custom Procedure packs.

Pennine's vision is to be a "world class, innovative manufacturer, supplying the global healthcare market with outstanding commitment to its people, customers and suppliers".

Pennine is enormously proud of its achievements but knows that it is only successful because of its people and real values, which are at the heart of everything they do

To be progressive
To be responsive
To be passionate
To have integrity

With a culture and history of innovation and engineering expertise, Pennine is still a competitive supplier of disposables on a worldwide basis and today embraces the skills and techniques of 'Lean Sigma' manufacture, ensuring that efficiency and waste reduction are targeted and that Pennine staff are highly trained, motivated and entirely focussed on the customer.

Pennine is consistently pursuing new opportunities and embracing innovative solutions, working with clinicians and healthcare professionals to develop new products for advanced patient care.

For further information, please call our sales team on +44 (0) 1332 794880 or email us at: marketing@penninehealthcare.co.uk. Alternatively, you can view all our product ranges at our website www.penninehealthcare.co.uk



Microsuction Cannulae



Urology



Custom Procedure Packs



Respiratory Products



Surgical Suction



Gastroenterology

Pennine[®]
H E A L T H C A R E

Pennine Healthcare

City Gate
London Road
Derby DE24 8WY
United Kingdom
Telephone: +44 (0) 1332 794880
Fax: +44 (0) 1332 794890

Email: info@penninehealthcare.co.uk
Web: www.penninehealthcare.co.uk

Some of the products in this catalogue are subject to a minimum order quantity.
Please ask our sales department for further details.

Ivor Shaw Ltd trading as Pennine Healthcare