




Ligating Clips

Ligating Clips



Endoscopic Appliers

CE 0473



Ligating Clips

Reference Code	Size	Compatible with applicators (box jaws grip mm)	Product Specification	Colour	Pcs per cartridge	Cartridges per box
WM-100	Small	0.59-0.65	V-shaped titanium ligating clips	Blue	6	36
WM-200	Medium	0.84-0.90		White	6	36
WM-300	Medium-Large	1.20-1.30		Green	6	20
WM-400	Large	1.30-1.50		Yellow	6	20

Open Appliers

Reference Code	Clip size	Applier length	Pcs per box
WM-S18	Small	18	1
WM-S28	Small	28	1
WM-M23	Medium	23	1
WM-M28	Medium	28	1
WM-ML28	Medium-Large	28	1
WM-L28	Large	28	1

Endoscopic Appliers - Front loaded clip applicator

Reference Code	Clip size	Pcs per box
WM-205ML	Medium-Large	1
WM-205L	Large	1

Endoscopic Appliers - Side loaded clip applicator

WM-305ML	Medium-Large	1
----------	--------------	---