



DISPOSABLE HEALTHCARE PRODUCTS



Microsuction Cannulae

WORKING IN PARTNERSHIP WITH
HEALTHCARE PROFESSIONALS

Pennine[®]
HEALTHCARE



Contents

Page

Zoellner Range

Zoellner Suction Probes	3
Zoellner Soft Tip	3
Zoellner with Vacuum Control	3
Zoellner with Flexiset	3

Fine Ends and Bellucci Range

Fine Ends	4
Bellucci Suction Probes	4
Bellucci with Flexiset	4

Frazier Range


Frazier Suction Probes	5
Frazier Soft Tip	5

Laryngeal and Oesophageal Tubes

Laryngeal Suction Probes	6
Oesophageal Suction Probes	6

Sundry

Sundry Suction Probes	7
Adson Suction Probes	7
Magill Suction Probes	7
Fess Suction Probes	7
Schuknecht Suction Probes	7
Barron Suction Probes	7
Flexiset	7
Adaptor	7



Zoellner Suction Probes

Single-use suction probes designed for ENT and any surgical procedure where fine suction is required.

- Sterile, individually wrapped, peel pouch
- Stainless steel probe
- Ergonomic handle
- 10° and 30° angle option
- Distal bend

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NZ-0001	Plain Tip	85mm	30°	50	6FG	2.0mm
NZ-0002	Plain Tip	85mm	10°	50	6FG	2.0mm

- Integral atraumatic soft tip to reduce trauma in delicate areas

Zoellner Soft Tip

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NZ-0001/S	Soft Tip	85mm	30°	50	6FG	2.0mm
NZ-0002/S	Soft Tip	85mm	10°	50	6FG	2.0mm

- Fingertip control for precision suction


Zoellner with Vacuum Control

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NZ-0008	Plain Tip	85mm	30°	50	6FG	2.0mm

- Towel clip for drape attachment
- Step adaptor for main suction
- Reduces fatigue due to tubing drag

Zoellner with Flexiset

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NZS-0002	Plain Tip	85mm	10°	50	6FG	2.0mm





Fine Ends

Precision bore stainless steel Fine Ends for use with Zoellner and Bellucci suckers.

- Sterile, individually wrapped, peel pouch
- Colour coded for easy identification
- Precision bore

Code No.	Type	Length	Qty. Per Box	Gauge	Ømm
NZ-0003	Fine End	40mm	50	18 gauge	1.27mm
NZ-0004	Fine End	40mm	50	20 gauge	0.91mm
NZ-0005	Fine End	40mm	50	22 gauge	0.72mm
NZ-0006	Fine End	40mm	50	24 gauge	0.57mm
NZ-0007	Fine End	40mm	50	26 gauge	0.46mm

Bellucci Suction Probes

Single-use suction probes suitable for use in theatres and in high risk procedures. Vacuum controlled for precise suction control and 'Sure-Grip' handle for comfort and secure handling.

- Sterile, individually wrapped, peel pouch
- Stainless steel probe
- Ergonomic 'Sure-Grip' handle
- Proximal bend with 30° angle
- With vacuum control

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NB-0015	Plain Tip	80mm	30°	50	6FG	2.0mm




Bellucci with Flexiset

Bellucci suction probe pre-attached to flexible, lightweight tubing.

- Towel clip for drape attachment
- Step adaptor for main suction
- Lightweight tubing reduces fatigue caused by tubing drag
- Proximal bend with 30° angle
- With vacuum control

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NBS-0015	Plain Tip	80mm	30°	50	6FG	2.0mm





Frazier Suction Probes

Single-use probes manufactured with stainless steel and 'Sure-Grip' handle for ease of use.


- Particularly suitable for neuro and spinal surgery where delicate manipulation is required.
- Sterile, individually wrapped, peel pouch
- Vacuum controlled for precise suction control
- Includes flexible stylette for tube clearance

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NF-0009	Plain Tip	120mm	30°	50	6FG	2.0mm
NF-0011	Plain Tip	120mm	30°	50	9FG	3.0mm
NF-0013	Plain Tip	120mm	30°	50	12FG	4.0mm

- Available with atraumatic soft tip for use in delicate areas

Frazier Soft Tip

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NF-0010	Soft Tip	120mm	30°	50	6FG	2.0mm
NF-0012	Soft Tip	120mm	30°	50	9FG	3.0mm
NF-0014	Soft Tip	120mm	30°	50	12FG	4.0mm





Laryngeal Suction Probes

Suitable for laryngoscopy, pharyngoscopy and bronchoscopy procedures. Varying lengths and gauges are available depending on the procedure.


- Sterile, individually wrapped, peel pouch
- Vacuum controlled for precise suction control
- Atraumatic soft tip

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NL-0023	Soft Tip	250mm	80°	10	6FG	2.0mm
NL-0024	Soft Tip	350mm	80°	10	6FG	2.0mm
NL-0025	Soft Tip	250mm	80°	10	9FG	3.0mm
NL-0026	Soft Tip	350mm	80°	10	9FG	3.0mm

Oesophageal Suction Probes

Code No.	Type	Length	Angle	Qty. Per Box	CH/FG	Ømm
NO-0027	Soft Tip Oesophageal	540mm	80°	10	12FG	4.0mm





Sundry Suction Probes

Full range of supplementary suction probes for a variety of different procedures based on individual requirements.

- Magill and Adson suction probes are suitable for maxfax and dental procedures
- Fess suction probes are specific to functional endoscopic sinus surgery
- Schuknecht and Barron suction probes are suitable for ear surgery procedures where delicate manipulation is required
- Sterile, individually wrapped, peel pouch
- Vacuum controlled for precise suction control

Code No.	Type	Length	Angle	Qty. Per Box	Size	Ømm
NA-0016	Adson Suction Probes	120mm	30°	50	9FG	3.0mm
NT-0017	Magill Suction Probes	120mm	5°	50	9FG	3.0mm
NF-0018	Fess Suction Probes, Soft Tip	120mm	80°	50	18G	1.27mm
NF-0019	Fess Suction Probes, Soft Tip	120mm	80°	50	16G	1.65mm
NS-0020	Schuknecht Suction Probes	90mm	30°	50	20G	0.91mm
NS-0021	Schuknecht Suction Probes	90mm	30°	50	24G	0.57mm
NBA-0001	Barron Suction Probes	80mm	30°	50	18G	1.27mm
NBA-0002	Barron Suction Probes	80mm	30°	50	16G	1.65mm

Flexiset

Highly flexible lightweight tubing with a unique towel clip which secures the tubing to the drape.

Code No.	Type	Length	Qty. Per Box
NSS-0022	Flexiset with Towel Clip	620mm	50



Adaptor

Enables non-vacuum sucker to be used with vacuum control.

Code No.	Type	Length	Qty. Per Box
NS-0022	House Adaptor	55mm	50

Please note: All angles and lengths are nominal

Tel: 44 (0) 1332 794880 Email: marketing@penninehealthcare.co.uk www.penninehealthcare.co.uk

LATEX
FREE



Pennine Healthcare is one of the UK's largest privately owned manufacturers of single-use sterile medical devices, with a customer base, not only in the UK but in over 50 countries worldwide.

Pennine has an extensive range of disposables for use in urology, anaesthetics, respiratory, gastroenterology and surgical suction, as well as being a leading supplier of Custom Procedure packs.

Pennine's vision is to be a "world class, innovative manufacturer, supplying the global healthcare market with outstanding commitment to its people, customers and suppliers".

Pennine is enormously proud of its achievements but knows that it is only successful because of its people and real values, which are at the heart of everything they do

To be progressive
To be responsive
To be passionate
To have integrity

With a culture and history of innovation and engineering expertise, Pennine is still a competitive supplier of disposables on a worldwide basis and today embraces the skills and techniques of 'Lean Sigma' manufacture, ensuring that efficiency and waste reduction are targeted and that Pennine staff are highly trained, motivated and entirely focussed on the customer.

Pennine is consistently pursuing new opportunities and embracing innovative solutions, working with clinicians and healthcare professionals to develop new products for advanced patient care.

For further information, please call our sales team on +44 (0) 1332 794880 or email us at: marketing@penninehealthcare.co.uk. Alternatively, you can view all our product ranges at our website www.penninehealthcare.co.uk



Gastroenterology



Urology



Custom Procedure Packs



Respiratory Products



Surgical Suction



Airways Management

Pennine[®]
H E A L T H C A R E

Pennine Healthcare

City Gate
London Road
Derby DE24 8WY
United Kingdom
Telephone: +44 (0) 1332 794880
Fax: +44 (0) 1332 794890

Email: info@penninehealthcare.co.uk
Web: www.penninehealthcare.co.uk

Some of the products in this catalogue are subject to a minimum order quantity.
Please ask our sales department for further details.

Ivor Shaw Ltd trading as Pennine Healthcare