




# Reloadable Linear Cutter Stapler



## Disposable Stapler with Titanium Staples

CE 0473

Instrument Code	Cartridge Code	Size	Tissue Thickness	Cartridge Colour
WM-LC-60-38	WM-LCC-60-38	60mm	1.5	Blue
WM-LC-60-48	WM-LCC-60-48	60mm	2.0	Green
WM-LC-80-38	WM-LCC-80-38	80mm	1.5	Blue
WM-LC-80-48	WM-LCC-80-48	80mm	2.0	Green
WM-LC-100-38	WM-LCC-100-38	100mm	1.5	Blue
WM-LC-100-48	WM-LCC-100-48	100mm	2.0	Green