



## Insufflation Tubing

For use in laparoscopic surgery (minimally invasive) for delivery of clean, filtered CO<sup>2</sup>

This system is to be used only for the insufflation of carbon dioxide into the abdomen during laparoscopy.

- 0.2µm hydrophobic filters prevents cross contamination from patient to machine and vice versa
- 6mm diameter lightweight tubing in 2.0, 2.5 and 3.0 metre lengths
- Choice of machine end connectors for compatibility with standard insufflators
- Sterile, double-wrapped individual peel pouches

CE 0473

ISO 9001:2008 certification  
ISO 13485:2003 certification

**LATEX  
FREE**



INVESTOR IN PEOPLE

|   | Code No.    | Fittings   | Tubing Diameter | Length of Tubing |
|---|-------------|--|-----------------|------------------|
|    | INS-1320/LL | Male Luer / Luer Lock Connector<br>Fit to Wolf, Storz and Linvatec insufflators.     | 6mm             | 2m               |
|    | INS-1330/LL | Male Luer / Luer Lock Connector  | 6mm             | 3m               |
|    | INS-1330/E  | Luer Lock / Funnel Fit Connector<br>Fits to Olympus, Storz and Stryker insufflators. | 6mm             | 3m               |
|  | INS-1250    | Male Luer Lock / Wolf Connector<br>Fits to Wolf, Storz and Linvatec insufflators.    | 6mm             | 2.5m             |
|  | INS-1332    | Luer Lock / Push Fit Connector   | 6mm             | 3m               |

Fits a wide range of insufflators such as:-

Wolf, Storz, Linvatec, Olympus, Stryker, Smith & Nephew and Snowden Pencer\*.

\*Please check compatibility with the appropriate manufacturer.

**Pennine Healthcare**

City Gate, London Road, Derby, DE24 8WY

Tel: +44 (0) 1332 794880

Fax: +44 (0) 1332 794890

