



Insufflation Tubing

For use in laparoscopic surgery (minimally invasive) for delivery of clean, filtered CO²

This system is to be used only for the insufflation of carbon dioxide into the abdomen during laparoscopy.

- 0.2µm hydrophobic filters prevents cross contamination from patient to machine and vice versa
- 6mm diameter lightweight tubing in 2.0, 2.5 and 3.0 metre lengths
- Choice of machine end connectors for compatibility with standard insufflators
- Sterile, double-wrapped individual peel pouches

CE 0473

ISO 9001:2008 certification
ISO 13485:2003 certification

**LATEX
FREE**



INVESTOR IN PEOPLE

	Code No.	Fittings	Tubing Diameter	Length of Tubing
	INS-1320/LL	Male Luer / Luer Lock Connector Fit to Wolf, Storz and Linvatec insufflators.	6mm	2m
	INS-1330/LL	Male Luer / Luer Lock Connector	6mm	3m
	INS-1330/E	Luer Lock / Funnel Fit Connector Fits to Olympus, Storz and Stryker insufflators.	6mm	3m
	INS-1250	Male Luer Lock / Wolf Connector Fits to Wolf, Storz and Linvatec insufflators.	6mm	2.5m
	INS-1332	Luer Lock / Push Fit Connector	6mm	3m

Fits a wide range of insufflators such as:-

Wolf, Storz, Linvatec, Olympus, Stryker, Smith & Nephew and Snowden Pencer*.

*Please check compatibility with the appropriate manufacturer.

Pennine Healthcare

City Gate, London Road, Derby, DE24 8WY

Tel: +44 (0) 1332 794880

Fax: +44 (0) 1332 794890

